



REHABILITATION TREADMILL

REHAB 403



Oceanic Rehabilitation Treadmill represents the top in medical rehabilitation technology. Its focus is to increase the strength of lower extremity walking speed and endurance. Combining forward and backward walking improves postural control and gait performance and lowers cardiovascular disease risk factors for stroke survivors and also prevention of future cardiovascular and cardiovascular events.

QUICK FACTS

Ac Motor	Model Rehab 403 - 4 HP AC
Set-Up Dimensions	L 82" x W 34" x H x 58"
Belt Size	3500 mm x 475 mm
Starting Increment 0.1 km/hr	Yes
Incline	Up to 15%
Decline Option	Yes
LED Screen	Yes
Reverse walking	Yes
Handle (with height and width adjustable side rail)	Yes
Frame Section	4" x 2"
Driver Rollers	114 mm
Power	220 Volts A.C, 20 Amps, 50Hz
Running Surface	63" x 18.5"
Product Weight	250 Kg
User Weight	180 kg
Therapist can Stand with Patient	Yes
Ramp	Yes

Oceanic Rehab treadmill is for slow walking with speed as slow as .01 km/hr and increment as low as 0.1 km/hr. We suggest using this Rehab treadmill with our Unweighing System OCU02 with Harness support for enhanced safety. Used to restore mobility in patients suffering from Stroke, Paralysis or fractures it will help to restore muscle strength, physical performance and for gait training. Used under the guidance of physical therapists and in Physical therapy centers.



OCEANIC FITNESS

B-36 INDUSTRIAL AREA
CHANALON, KURALI
PUNJAB
(An ISO 9001-2008 certified company)

WWW.OCEANICFITNESS.COM
Email: oceanicfitness@gmail.com
Mobile/WhatsApp +91 9872414100, 9814148008